

FOR IMMEDIATE RELEASE:

CONTACT:

Bridget Harrington

KUBTEC

Tel: (203) 364 – 8544

Fax: (203) 255 – 7494

Email: kubtec@kubtec.com

Website: www.kubtec.com

KUBTEC Awarded FDA Approval of DIGIVIEW 250 For Neonatal Imaging

Milford, CT, October 25, 2011 -- KUBTEC, the leader in innovative technology for digital x-ray systems, is pleased to announce FDA approval of their DIGIVIEW 250®, the first and only digital radiography detector that employs CMOS technology. This breakthrough technology provides low-dose and high-resolution imaging. The DIGIVIEW 250 CMOS digital radiographic detector is ideal for use in neonatal, pediatric and extremity imaging.

Weighing less than 7 lbs. and configured in a durable, custom-molded plastic housing with a integrated handle for easy positioning, the DIGIVIEW 250 CMOS digital radiographic detector features 96µm pixel and 14 bit contrast with a 20 cm x 25 cm active imaging area, providing the most detailed, high-resolution and high-contrast images. For added ease of use, the DIGIVIEW 250 integrates with any existing x-ray source instantly without hardware connections. Designed for the neonatal environment, but suitable for a wide-range of low-dose applications, the DIGIVIEW 250 acquires detailed images with lower radiation dose. For neonatal applications, the DIGIVIEW 250's compact format and easy-to-position design speeds imaging while greatly reducing radiation exposure.

With rapid image acquisition and display, KUBTEC's DIGIVIEW 250 with its full complement of DICOM 3.0 compliant, user-friendly image enhancement tools, on screen measurement capabilities, ability to annotate and send images to all PACS systems, the DIGIVIEW 250 brings high speed, high resolution digital imaging to the most fragile patients.

To learn more about the DIGIVIEW 250 CMOS digital radiographic detector, please visit KUBTEC's booth (#5517) at the RSNA's (Radiological Society of North America) annual meeting in Chicago, November 27 – December 1, 2011. Our staff will be happy to discuss the advanced features of this breakthrough system with you.

#####

ABOUT KUBTEC – KUBTEC and DIGIVIEW are registered trademarks of KUB Technologies, Inc. a manufacturer of compact digital x-ray inspection systems for a variety of applications including specimen radiography, science and research, irradiation, security and non-destructive testing. Systems and components along with local support are available around the world.