

---

## FOR IMMEDIATE RELEASE:

### CONTACT:

**Bridget Harrington**

**KUBTEC**

**Tel: (203) 364 – 8544**

**Fax: (203) 255 – 7494**

**Email: [kubtec@kubtec.com](mailto:kubtec@kubtec.com)**

**Website: [www.kubtec.com](http://www.kubtec.com)**

## **Texas Children’s Hospital Purchases KUBTEC’s XPERT 80-L Digital Specimen Radiography System for Anatomic Pathology**

Milford, CT, October 7, 2011-- KUBTEC, the innovator and manufacturer of digital x-ray technology, announces the sale of their XPERT 80-L unit with a 17” x 17” detector, to Texas Children’s Hospital, Houston. Texas Children’s Hospital is the largest children’s hospital in the nation, and in 2011 it was named as one of the top four Pediatric Facilities in the United States by U.S. News and World Report. Texas Children’s Hospital is also affiliated with the Baylor College of Medicine, as their main pediatric training site.

Through their affiliation with the Baylor College of Medicine, the decision makers at Texas Children’s were able to see the image quality, the speed and ease of operation of another KUBTEC system that had been installed and in use for a few years at Baylor. With the largest detector available in a cabinet x-ray system, the superior image quality, rapid acquisition and user-friendly image enhancement software, and the endorsement of the users at Baylor College of Medicine, the decision was clear. Customer support and service were also big factors in selecting the KUBTEC XPERT 80-L, said Felton Nails, Histology Manager. Known for their responsive customer support and servicing, KUBTEC was the obvious choice.

A driving force in innovation for cabinet x-ray technology, KUBTEC’s XPERT family of digital x-ray systems is available with digital detectors from 2” x 2” up to 17” x 17” and with optional x-ray sources. With up to 16 bit contrast resolution KUBTEC systems acquire the sharpest images and have the ability to penetrate the densest specimens. All the XPERT systems come with a number of software features; AEC, auto-calibration, automatic window leveling, and a full suite of user-friendly image enhancements, including 3D, measurement tools and the ability to save images in multiple formats. KUBTEC systems also provide touch screen technology and interface seamlessly with most Radiology, Pathology and Research management software suites.

#####

**ABOUT Kubtec** – Kubtec is a registered trademark of KUB Technologies, Inc. a manufacturer of compact digital x-ray inspection systems for a variety of applications including specimen radiography, science and research, irradiation, security and non-destructive testing. Systems and components along with local support are available around the world.