
FOR IMMEDIATE RELEASE:

CONTACT:

Cindy Cavoto

KUBTEC

Tel: (203) 364 – 8544

Fax: (203) 255 – 7494

Email: kubtec@kubtec.com

Website: www.kubtec.com

KUBTEC introduces new Specimen Radiography System for the Biopsy Suite at RSNA 2008

Milford, CT, September 26, 2008. Kubtec will introduce their newest product innovation, the XPERT 20 at RSNA 2008. A digital specimen radiography system exclusively for the biopsy suite, the XPERT 20 is available with either a 50mm x 50mm detector or a 100mm x 50mm detector in a self-contained, shielded cabinet x-ray system that takes up less space than a PC monitor.

A small focal spot x-ray source provides increased visibility and sharpness of detail for lesions and abnormalities in breast tissue samples. With an X-ray source of up to 25kV and 1.0mA, the XPERT 20™ provides images with the highest resolution to locate the smallest micro-calcifications.

All systems include DICOM compliant, user-friendly software with AEC (automatic exposure control), auto-calibration, automatic window leveling, the ability to send multiple images with annotations to PACS, and zero warm up time. Where space is at a premium and efficient, high quality patient care is the goal, the low cost digital XPERT 20™ is the economical solution.

#####

ABOUT Kubtec – Kubtec is a registered trademark of KUB Technologies, Inc. a manufacturer of compact digital x-ray inspection systems for a variety of applications including specimen radiography, science and research, irradiation, security and non-destructive testing. Systems and components along with local support are available around the world.