
FOR IMMEDIATE RELEASE:**CONTACT:****Bridget Harrington****KUBTEC****Tel: (203) 364 – 8544****Fax: (203) 255 – 7494****Email: kubtec@kubtec.com****Website: www.kubtec.com**

HHMI/Baylor College Medicine takes delivery of XPERT 80

Milford, CT, August 24, 2011 -- KUBTEC, the innovator in digital cabinet x-ray systems designed for clinical and research applications, announce the delivery of their XPERT 80 to HHMI/Baylor College of Medicine. The XPERT 80 will support research in the laboratory of Brendan Lee, M.D., Ph.D, a professor in the Department of Molecular and Human Genetics. Dr. Lee also holds the distinction of being named a Howard Hughes Medical Institute (**HHMI**) Investigator. In addition to being a professor in Baylor's Department of Molecular and Human Genetics and an HHMI Investigator, he is the Chief of the Skeletal Dysplasia Clinic at Texas Children's Hospital and the director of the Bone Disease Program of Texas.

Jianning Tao, Ph.D, Postdoctoral Fellow and Principal Researcher, in Dr. Lee's HHMI laboratory at Baylor College of Medicine, said, in the justification to purchase, "The XPERT 80 has been chosen after a thorough investigation of the available products on the market. It is a top-tier specimen radiography system and is the only system available with 16-bit contrast resolution and 1.0 mA x-ray source enabling it to image the tiniest micro-calcifications in the specimens."

A driving force in innovation for cabinet x-ray technology, KUBTEC's XPERT family of digital x-ray systems is available with digital detectors from 2" x 2" up to 17" x 17" and with optional x-ray sources up to 160kV. With up to 16 bit contrast resolution KUBTEC systems acquire the sharpest images and have the ability to penetrate the densest samples. All the XPERT systems come with a number of software features; AEC, auto-calibration, automatic window leveling, and a full suite of user-friendly image enhancements, including 3D, measurement tools and the ability to save images in multiple formats. KUBTEC systems also provide up to 5X imaging capabilities and touch screen technology.

#####

ABOUT Kubtec – Kubtec is a registered trademark of KUB Technologies, Inc. a manufacturer of compact digital x-ray inspection systems for a variety of applications including specimen radiography, science and research, irradiation, security and non-destructive testing. Systems and components along with local support are available around the world.